

Adult SCRAM™ Bag Kit List



Find out more about the features of Adult SCRAM™ at
www.scramsystem.org

Use this List as a Guide only.
Review your requirements and add or delete items to meet your needs.

Adult SCRAM™ - Example Kit List

Primary Area	Quantity	Checked
MacIntosh blade size 4	1	
MacIntosh blade size 3	1	
Adult Laryngoscope Handle	2	
20ml syringe	1	
Tube tie - check length 100cm	2	
Lubricating gel	2	
Adult Catheter Mount (In packaging)	1	
Adult Filter (In packaging)	1	
End-Tidal CO2 (ETCO2) sampling line	1	
Adult Magill Forceps	1	
ET Tube 8.0 cuffed uncut	1	
ET Tube 7.0 cuffed uncut	1	
Large Stylet	1	
Medium Stylet	1	
Bougie 10Ch	1	
Bougie 15Ch	1	
Teams who do not mind having a gradual bend in the bougie store them in the primary use area along the spine of the bag, whereas teams that prefer a straight bougie store them in a separate bougie holder.		
Secondary Area	Quantity	Checked
MacIntosh blade size 2	1	
Guedel Airways 4,3,2	1 each	
Guedel Airways 1,0	1 each	
Nasal Airways 7.0 & 6.0	1 each	
Video Laryngoscope blade size 3	1	
Video Laryngoscope blade size 4	1	
Tape zinc oxide	1 roll	
PEEP Valve	1	

Surgical Airway Pouch	Quantity	Checked
Tracheostomy Kit	1	
Tracheal Hook	1	
Checklist Pouch	Quantity	Checked
RSI Checklist & Rapid Induction Checklist	1 each	
Skin Marker Pen (<i>in loop below checklist pouch</i>)	1	
File Storage System	Quantity	Checked
Orogastric tube 16F	1	
Drainage bag	1	
OG Tube Spigot	1	
Second Generation Supraglottic Airway Device	1	
Nasal Cannula	1	
Orange Disp Bag	1	
Pouch with Yellow Tag	Quantity	Checked
Scalpel Size 10	1	
ET Tube 6.0 cuffed uncut	1	
Pouch with Blue Tag (Thoracostomy kit)	Quantity	Checked
Scalpel size 24	1 each	
Spencer Wells 20cm	1	
Chloroprep 3ml	1	

Personalised kit

Personalised kit		



The Example Kit List for the **Adult SCRAM™ bag** is intended only as an example of how the Adult SCRAM™ bag may be stocked. The Adult SCRAM™ bag may be stocked in different ways to accommodate the requirements of different users. You should review your requirements carefully and add or delete items to meet those needs. It is the responsibility of the health professional or health organisation stocking the SCRAM™ bag to ensure that any equipment and other items stocked in it meet with and comply with its own requirements and all and any relevant guidelines, statements of practice or laws in the location where the SCRAM™ bag will be used. The Example Kit List does not constitute advice and no representations or warranties, express or implied, are given that the Example Kit List is accurate or complete or in line with current medical guidelines or practice. No liability will accrue to the manufacturer, distributor, **Scottish Health Innovations Limited (SHIL)** or any other party in the event that the User suffers any loss or damage as a result of relying upon the suggestions made in the Example Kit List, or any omissions from it.